



Serum Triglyceride Assay Kit

Cat# STG-1-NC

INSTRUCTION MANUAL ZBM0036.06

STORAGE CONDITIONS

- Diluent & Buffers
 - Room Temperature
- Glycerol Reagent A:
 - +4°C **Use reconstituted Glycerol Reagent A within 7 days.**
- Glycerol Standard & Reagent B
 - -20°C **NEW REAGENT B STORAGE TEMPERATURE**

ALL ZEN-BIO INC PRODUCTS ARE FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE OR FOR USE IN DIAGNOSTIC OR CLINICAL PROCEDURES.

ORDERING INFORMATION AND TECHNICAL SERVICES

- Zen-Bio, Inc.
- 3920 South Alston Avenue
- Durham, NC 27713, USA
- Telephone (919) 547-0692
- Facsimile (FAX) (919) 547-0693
- Toll Free 1-866-ADIPOSE (866)-234-7673
- Electronic mail (e-mail) information@zen-bio.com
- World Wide Web <http://www.zenbio.com>

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, expressed or implied, including without limitation implied warranties of merchantability or fitness for a particular purpose, are provided by Zen-Bio, Inc. Zen-Bio, Inc. shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.




Contact ZenBio, Inc. within no more than 24 hours after receipt of products for all claims regarding shipment damage, incorrect ordering or other delivery issues. Delivery claims received after 7 days of receipt of products are not subject to replacement or refund.

INTRODUCTION

This kit is designed to accurately determine the amount of triglyceride present in blood serum or plasma of humans, mice, rats, and other animals in a 96-well format for increased throughput analysis. Blood can be collected in plain evacuated tubes or in the presence of common anti-coagulants: heparin, EDTA, sodium citrate, or ammonium oxalate. Serum should be separated from clotted blood by centrifugation as soon as possible and may be stored frozen (-20°C) prior to analysis.

PRINCIPLE OF THE ASSAY

Enzymatic hydrolysis of serum triglycerides by lipase generates free fatty acids and glycerol. The glycerol released is phosphorylated by adenosine triphosphate (ATP) forming glycerol-1-phosphate (G-1-P) and adenosine-5'-diphosphate (ADP) in the reaction catalyzed by glycerol kinase. G-1-P is then oxidized by glycerol phosphate oxidase to dihydroxyacetone phosphate (DAP) and hydrogen peroxide (H₂O₂). A quinoneimine dye is produced by the peroxidase catalyzed coupling of 4-aminoantipyrine (4-AAP) and sodium N-ethyl-N-(3-sulfopropyl)m-anisidine (ESPA) with H₂O₂, which shows an absorbance maximum at 540nm. The increase in absorbance at 540nm is directly proportional to glycerol (and triglyceride) concentration of the sample.



ITEMS INCLUDED IN THE KIT

The contents of this kit are sufficient for the assay of 1, 96 well plate. For additional reagents or formats, please order our bulk kit (cat# STG-5RB).

ITEM	DESCRIPTION	UNIT	QTY	STORAGE
Reagent A	Reconstitute w/ 11.0 ml deionized water prior to use. Use reconstituted reagent within 7 days!	BOTTLE	1	4°C
Reagent B	2.5 ml Reagent B solution new storage temperature	BOTTLE	1	-20°C
Glycerol standard	10 µl Glycerol @ 10mM [Reconstitute with 40 µl Standards Diluent to make the 2 mM glycerol standard; see page 4 for recommended dilution scheme]	VIAL	1	-20°C
Diluent	Standards Diluent, 2.0ml	BOTTLE	1	4°C
Tray	Clear polyvinyl tray for use with multi-channel pipetters	EACH	2	ROOM TEMP.
Assay Plates, blank	96-well assay plate, blank	PLATE	2	ROOM TEMP.


Other equipment/reagents required but not provided with the kit:

- ▶ Single-channel pipetter
- ▶ Multi-channel pipetter
- ▶ Plate reader with a filter of 540 nm
- ▶ Tubes for diluting glycerol standards

REAGENT PREPARATION

1. Thaw Reagent B to room temperature and gently invert. DO NOT VORTEX!
Reagent B may be refrozen (-20°C) and stored. Avoid repeated freeze/thaws.
2. Prepare the Glycerol Reagent A by adding 11.0 ml deionized water per bottle and gently invert. DO NOT VORTEX! Use a pipet to insure that the powder is completely dissolved. Keep at room temperature. If using a Reagent A solution previously prepared and stored at 2-8°C, also bring to room temperature. Store in a light protected bottle. Reconstituted Glycerol Reagent A is stable for 7 days refrigerated (2-8°C); store any remaining solution refrigerated (2-8°C).

3. Prepare Glycerol Standards using the 10mM glycerol standard provided. Briefly centrifuge thawed 10mM glycerol standard to collect all liquid in the bottom of the tube. Add 40 μl of diluent to 10mM stock glycerol standard to generate a 2mM glycerol standard. Add 25 μl of diluent to 6 tubes (not provided). Prepare the dilution series as outlined below. Mix each new dilution thoroughly before proceeding to the next. The 2 mM solution serves as the highest standard, and the diluent serves as the zero standard.



ASSAY PROCEDURE

STEPWISE GLYCEROL AND TRIGLYCERIDE MEASUREMENT

1. Warm reconstituted Glycerol Reagent A and Reagent B to room temperature.
2. Add 80 μl of reconstituted Glycerol Reagent A to the appropriate number of wells for samples and standard curve in the provided 96-well plate.
3. Add 5.0 μl of prepared standard curve dilutions to wells containing Glycerol Reagent A. Transfer 1.0 - 5.0 μl of serum or plasma samples to remaining wells. Incubate at room temperature for 15 minutes
4. Read at 540 nm using a microtiter plate reader and record absorbance values as Glycerol Readings. (Glycerol standards from this reading are used to calculate total glycerol)
5. Add 20 μl of Reagent B to each well and incubate and additional 15 minutes at room temperature. Record absorbance values as Total Triglyceride Readings. (Glycerol standards from this reading are used to calculate combined glycerol and triglyceride levels.)

DATA ANALYSIS


STEPWISE GLYCEROL AND TRIGLYCERIDE MEASUREMENT

Generate standard curves: see example below

[DO NOT use this standard curve to generate your data. This is an example.]

Subtract the OD value of the 0µM standard from all OD values including the standard curve.

µM Glycerol	Gly Std OD	TG std OD	Gly - Blank	TG -Blank
2000	0.293	0.342	0.254	0.299
1000	0.151	0.178	0.112	0.135
500	0.091	0.109	0.052	0.066
250	0.062	0.072	0.023	0.029
125	0.049	0.057	0.01	0.014
62.5	0.041	0.048	0.002	0.005
31.25	0.039	0.044	0	0.001
0	0.039	0.043		



	Gly	TG
slope =	0.0001	0.0001
intercept=	0.0065	0.0055
r²=	0.9965	0.9978

y = observed O.D. minus the blank

x = concentration of glycerol in µM

To calculate x for each y, (i.e. to change the observed O.D. into glycerol concentration) use the following equation:

y=(slope) times (x) plus intercept

y=mx+b so x=(y-b)/m

x=(y – 0.0065)/0.0001 where 0.0001= slope of the line and 0.0065= y intercept. Be careful to enter the proper sign for the y intercept value as it may be a negative number.

Any OD values greater than the highest standard (2000 µM) should be suspect. The compound should be re-assayed using a lower dose of the compound at treatment OR a dilute solution of the condition medium at the time of the assay.

The R² value should be equal or greater than 0.98 for the standard curve to be valid. Any R² values below 0.98 must have the standard curve run again.

Solve for the Total Glycerol concentration (i.e. total triglyceride concentration) for each OD. **Remember to include any Dilution Factor in the equation.** Data is expressed as µM Glycerol.

NOTE: Any OD values that are negative after the blank is subtracted should be considered to be 0 for the OD value.

1. Calculate the values for glycerol present in sample by using the initial glycerol standards reading.
2. Calculate the values for the combined glycerol and triglyceride by using the second glycerol standards reading.
3. Subtract the glycerol levels from step 1 from the combined glycerol and triglyceride values in step 2. This gives the value of total triglyceride in the sample.

TROUBLESHOOTING GUIDE

Problem	Suggestions
High background or the triglyceride reagent turns a darker color before the assay begins.	<ul style="list-style-type: none"> • Use clean tray and tips • Change pipet tips frequently
Edge effects	<ul style="list-style-type: none"> • Ensure a saturated humidity in the incubator to prevent evaporation from the outside wells
Inconsistent OD reading	<ul style="list-style-type: none"> • Be careful when pipetting to avoid bubbles. If bubbles persist, burst the bubbles using a large gauge needle prior to reading and read the plate again. • Mix the lysates well before transferring the 20μl to the Wash buffer plate.

REFERENCES

1. Green, H. and Kehinde, O. (1974) Sublines of mouse 3T3 cells that accumulate lipid. *Cell* 1, 113-116.
2. Hauner, H., et al., (1989) *J. Clin. Invest.*(84), 1663-1670.
3. Kuri-Harcuch W, Wise LS, Green H. (1978) Interruption of the adipose conversion of 3T3 cells by biotin deficiency: differentiation without triglyceride accumulation. *Cell* 14:53-58.

APPENDIX A: PLATE LAYOUT

H	G	F	E	D	C	B	A	
								1
								2
								3
								4
								5
								6
								7
								8
								9
								10
								11
								12

APPENDIX B: PROCEDURE FLOWCHART

Stepwise Glycerol and Triglyceride Detection

Reconstitute Glycerol Reagent A and thaw Reagent B according to instructions



Add 80 μ l/well reconstituted Glycerol Reagent A



Add 5 μ l/well glycerol standards and 1-5 μ l/well samples.




Incubate 15 minutes @ room temperature. Measure the optical density of each well at 540 nm using a spectrophotometer plate



Add 20 μ l/well Reagent B



Incubate 15 minutes @ room temperature. Measure the optical density of each well at 540 nm using a spectrophotometer plate reader.



Sample



$\xrightarrow{5 \mu\text{l}}$

