



Lipid Staining Kit

200 Point Assay Kit in 96-well Format

Cat# ST-R100

INSTRUCTION MANUAL ZBM0012.01

STORAGE CONDITIONS

- **Reagents & Buffers:** Room temperature 25°C

For *in vitro* Use Only

For Research Use Only. Not For Use In Diagnostic Procedures

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, expressed or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Zen-Bio, Inc. Zen-Bio, Inc. shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

- **Zen-Bio, Inc.**
- **3200 Chapel Hill-Nelson Blvd., Suite 104**
- **PO Box 13888**
- **Research Triangle Park, NC 27709**
- **Telephone** (919) 547-0692
- **Facsimile (FAX)** (919) 547-0693
- **Toll Free** 1-866-ADIPOSE (866)-234-7673
- **Electronic mail (e-mail)** information@zen-bio.com
- **World Wide Web** <http://www.zen-bio.com>

INTRODUCTION

This protocol is designed to stain cells in a 96-well format using Oil Red O. Oil Red O is a fat-soluble diazo dye that stains neutral triglycerides and lipids. For other formats, please adjust the volumes added to each well according to the surface area of the well/flask you are using.

ITEMS INCLUDED IN THE KIT

ITEM	STORAGE NOTES	UNIT	QTY
Fixative Solution	Store at room temperature, preferably in a cabinet approved for storage of flammable chemicals.	40 ML BOTTLE	1
Oil Red O Stock Solution	Store at room temperature in a dark bottle, preferably in a cabinet approved for storage of flammable chemicals.	6 ML BOTTLE	1
Syringe	10 cc syringe	EACH	1
Filter	1 μ m PTFE filter	EACH	1
Tray	For multi-channel pipettes, clear polyvinyl	EACH	3

Other equipment/reagents required but not provided with the kit:

- Mature adipocytes or other lipid-containing cells to be stained
- Deionized or distilled water (dH₂O)
- Multi-channel Pipet, single channel pipette and pipette tips
- Tubes for preparing working solutions

ASSAY PROCEDURE

96-well plate format Protocol

1. Remove row by row all the liquid and add 100 μ l/ well Fixative Solution. Incubate at room temperature for 5 minutes. Repeat liquid exchange with another 100 μ l/ well Fixative Solution.
2. Fix cells at least 1 hour at 4°C. Alternatively, the plate of cells in Fixative Solution may be sealed using sealing film/foil (not included) and stored at 4°C for 7 days before proceeding.
3. **Prepare Working Oil Red Solution (WOROS): (Prepared fresh on day of assay) 40% water / 60% Oil Red O stock.**
 - a. Add 6.0 ml of Oil Red O Stock Solution to 4.0 ml of water and **mix by inversion**.
 - b. This working solution must be kept at room temperature at least 20 minutes before filtering.
 - c. Filter out particulate material using a 1 μ m polypropylene PTFE syringe filter (supplied with the kit) or #1 Whatman filter paper (not provided). NOTE: Filter paper is effective; it just takes longer to filter. Polycarbonate filters will dissolve in isopropanol; DO NOT USE THEM.
4. Remove all the Fixative Solution. Wash cells at least twice with 200 μ l/well dH₂O using a multichannel pipette. Remove all the water using a multichannel pipette. Make sure all the water is removed **BEFORE** proceeding.
5. Add 50 μ l/ well Working Oil Red O solution (WOROS) at room temperature for 10-15 minutes. Be careful not to touch the sides of the wells; the pipette tips should go at a 90° angle into the wells.
6. Remove all of the Oil Red O solution from each well. Wash 2-3 times with 60-80 μ l/ well dH₂O. Pictures may be taken at this time. Leave water inside wells during photography.