

HUMAN ADULT DERMAL FIBROBLAST CARE MANUAL

INSTRUCTIONAL MANUAL ZBM0023.01

SHIPPING CONDITIONS

Human Adult Dermal Fibroblast Cells

Orders are delivered via Federal Express courier. All US and Canada orders are shipped via Federal Express Priority service and are usually received the next day. International orders are usually received in 3-4 days.

Must be processed upon shipment receipt.

STORAGE CONDITIONS

Media: Short Term 4°C 6 months -20°C

Cells: Frozen: liquid nitrogen Plated: 37°C incubator

All Zen-Bio Inc products are for research use only. Not approved for human or veterinary use or for use in diagnostic or clinical procedures.

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, expressed or implied, including without limitation implied warranties of merchantability or fitness for a particular purpose, are provided by Zen-Bio, Inc. Zen-Bio, Inc. shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Zen-Bio, Inc warrants its cells only if Zen-Bio media are used and the recommended protocols are followed. Cryopreserved human adult dermal fibroblasts are assured to be viable when thawed and maintained according to Zen-Bio protocols.

ORDERING INFORMATION AND TECHNICAL SERVICES

Zen-Bio, Inc.

3200 Chapel Hill-Nelson Blvd., Suite 104

PO Box 13888

Research Triangle Park, NC 27709

U.S.A.

Telephone

(919) 547-0692

Facsimile (FAX)

(919) 547-0693

Toll free (continental US only)

1-866-ADIPOSE 1-(866)-234-7673

Electronic mail (e-mail)

information@zen-bio.com

World Wide Web

<http://www.zen-bio.com>

HUMAN ADULT DERMAL FIBROBLASTS

<u>CONTENTS</u>	<u>PAGE #</u>
Introduction	3
Materials provided for each catalog item	3
Media composition	4
Pathogen testing	4
Plating and expansion of cryopreserved DF cells	5
Troubleshooting	6
Frequently asked questions	6

Introduction

The adult dermal fibroblast cells are isolated from the dermal layer of skin tissue from elective surgery in healthy non-diabetic donors between 18 and 60 years old. The cells are isolated either by centrifugal force following enzymatic treatment or from an explant culture. This instruction manual describes procedures to passage and culture the human dermal fibroblast cells.

Precautionary Notes

This product is for research use only. *It is not intended for human, veterinary, or in vitro diagnostic use.* Proper precautions and biological containment should be taken when handling cells of human origin, due to their potential biohazardous nature. **Always wear gloves and work behind a protective screen when handling primary human cells.** All media, supplements, and tissue cultureware used in this protocol should be sterile.

Human adult dermal fibroblast cell viability depends greatly on the use of suitable media, reagents, and sterile plastic wear. If these parameters are not carefully observed, limited differentiation may occur and cell growth may be slow.

Materials Provided for Each Catalog Item

- ❖ **Human Adult Dermal Fibroblasts** (96-,48-,24-,12-,6-well plates; 75-,25-cm² flasks)
 - Cat # DF-2096, -2048, -2024, -2012, -2006, -75, -25
- ❖ **Dermal Fibroblast Cryopreservation Medium**
 - Cat# DFM-100
- ❖ **Cryopreserved Human Adult Dermal Fibroblasts**
 - Cat # DF-F
 - Frozen vial containing 1 x10⁶ viable adult dermal fibroblasts (store in liquid nitrogen upon receipt)



Media Compositions

Dermal Fibroblast Culture Medium

cat # DF-1

500ml

DMEM
Fetal bovine serum
Penicillin
Streptomycin
Amphotericin B

Dermal Fibroblast Basal Medium

cat # DF-2

500ml

DMEM
Penicillin
Streptomycin
Amphotericin B

Dermal Fibroblast Freeze Medium

cat # DFM-100

100ml

DMEM
Fetal bovine serum
DMSO

NOTE: All dermal fibroblast media contains 4.15g/L D-glucose.

All Zen-Bio, Inc media are also available phenol red free.

Please inquire for custom media requests.


Pathogen Testing

Each donor is tested via PCR and found non-reactive for HIV-1, HIV-2, HTLV I, HTLV II, hepatitis B and hepatitis C. However, no known test can offer complete assurance that the cells are pathogen free. Our products are tested and are free from mycoplasma contamination. Proper precautions and biological containment should be taken when handling cells of human origin, due to their potential biohazardous nature. All human based products should be handled at a BSL-2 (Biosafety Level 2) or higher. Always wear gloves and work behind a protective screen when handling primary human cells.

Plating and Expansion Procedure for Cryopreserved Adult Dermal Fibroblasts

1. Remove cells from liquid nitrogen and place immediately into a 37°C water bath with agitation. Be careful not to submerge the cap of the vial into water. Do not leave the vials in water bath after most of the content has thawed. Rinse the vials with 70% ethanol before taking them to the culture hood.
2. Upon the thawing, add the cells to a sterile conical bottom centrifuge tube, containing 9 ml of Dermal Fibroblast Growth Medium (DF-1).
3. Centrifuge at 400 x g, 20°C, 10 minutes. Aspirate the medium and resuspend cells in a volume of DF-1 appropriate for counting the cells. Count using a hemacytometer.
4. Place approximately $0.75\text{--}1.0 \times 10^6$ cells in T-75 culture flasks using DF-1.
5. Incubate cells until they are 85-90% confluent (in about 4-5 days). Cells will need to be fed once a week with DF-1.
6. Aspirate medium and wash adult fibroblasts 4-5 times using sterile Phosphate Buffered Saline (PBS) to remove all traces of serum (until there is no foaming of the medium). Remove the PBS and release the cells from the flask bottom by adding 2 mL/T-75 flask (or 6 mL/T-225 flask) of 0.25% trypsin/ 2.21mM EDTA solution. Allow cells to trypsinize for 5 minutes at 37°C. Tap the flask gently to loosen the cells.
7. Neutralize the trypsin using 7 ml DF-1 per T-75 flask (or 21 ml per T-225 flask). Check the flask under a microscope to ensure all cells are free of the flask bottom.
8. Count the cells and plate in desired format. Ensure cells are evenly suspended when plating large numbers of plates or flasks. Do not agitate plates and flasks after plating. Place in a humidified incubator at 37°C and 5% CO₂, making sure the surface is level for even cell distribution.

ADULT DERMAL FIBROBLASTS



Troubleshooting Guide

Observation	Possible causes	Suggestions
Adult DF cells do not grow	<ol style="list-style-type: none"> 1. Cells have been passaged too many times 2. Cells expanded too high 	<ol style="list-style-type: none"> 1. Use cells of a lower passage number 2. Do not exceed 1:6 expansion ratio
Edge effects	<ol style="list-style-type: none"> 1. Medium in outside wells evaporated 	<ol style="list-style-type: none"> 1. Ensure a saturated humidity in the incubator and feed the cells no less than every 3 days. Make sure multiple plates are stacked no more than 3 plates high.

Frequently Asked Questions

- **Can I pass the cells?**
Dermal Fibroblast Cells can be trypsinized and replated several times. All cells are shipped after establishing a primary culture.
- **How fast do the cells replicate?**
The average doubling time is 18-24 hours. However, keep in mind that the replication rate for human dermal fibroblasts cells varies slightly from donor to donor.
- **Should antibiotics be included in the medium?**
Yes. Antibiotics and anti-fungal agents are always recommended since the cells are primary cells.
- **Where are the cells obtained?**
The adult dermal fibroblast cells are isolated from human dermal tissue.
- **Do you test for pathogens? Which ones?**
Yes. Each donor is tested via PCR and found non-reactive for HIV-1, HIV-2, HTLV I, HTLV II, hepatitis B and hepatitis C. However, since we cannot test all pathogens, please treat the culture as a potentially infectious agent.
- **What donor information do I receive?**
The donor's age, gender, and BMI are provided in the certificate of analysis that accompanies each lot of cells.
- **What is the formulation of Zen-Bio's serum-free media?**
Zen-Bio's serum-free media are not enhanced to supplement the absence of serum. These media are available for assay procedures where cells are rested from serum.