



# **96-well Serum/Plasma Glycerol Kit**

## **Free Glycerol Detection**

### **Cat# SGA-1**

**INSTRUCTION MANUAL    ZBM0043.00**

#### **STORAGE CONDITIONS**

---

- **Reagents & Buffers:** 4°C
- **Glycerol Standard:** -20°C
- **Blank assay plates (96-well):** Room Temperature

**For *in vitro* Use Only**

#### **LIMITED PRODUCT WARRANTY**

This warranty limits our liability to replacement of this product. No other warranties of any kind, expressed or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Zen-Bio, Inc. Zen-Bio, Inc. shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

#### **ORDERING INFORMATION AND TECHNICAL SERVICES**

- **Zen-Bio, Inc.**
- **3200 Chapel Hill-Nelson Blvd., Suite 104**
- **PO Box 13888**
- **Research Triangle Park, NC 27709**
- **Telephone** (919) 547-0692
- **Facsimile (FAX)** (919) 547-0693
- **Toll Free** 1-866-ADIPOSE (866)-234-7673
- **Electronic mail (e-mail)** [information@zen-bio.com](mailto:information@zen-bio.com)
- **World Wide Web** <http://www.zen-bio.com>

# INTRODUCTION

---



This kit is designed to accurately determine the amount of glycerol present in blood serum or plasma of humans, mice, rats, and other animals in a 96-well format for increased throughput analysis. Blood can be collected in plain evacuated tubes or in the presence of common anti-coagulants: sodium citrate, ammonium oxalate and EDTA. **NOTE: Heparin or Heparinized tubes should not be used because this will generate inaccurate readings.** Serum should be separated from clotted blood by centrifugation as soon as possible and may be stored frozen (-20°C) prior to analysis.

## PRINCIPLE OF THE ASSAY

---

### Detection of Free Glycerol

Glycerol present in sample is phosphorylated by adenosine triphosphate (ATP) forming glycerol-1-phosphate (G-1-P) and adenosine-5'-diphosphate (ADP) in the reaction catalyzed by glycerol kinase. G-1-P is then oxidized by glycerol phosphate oxidase to dihydroxyacetone phosphate (DAP) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). A quinoneimine dye is produced by the peroxidase catalyzed coupling of 4-aminoantipyrine (4-AAP) and sodium N-ethyl-N-(3-sulfopropyl)m-anisidine (ESPA) with H<sub>2</sub>O<sub>2</sub>, which shows an absorbance maximum at 540nm. The increase in absorbance at 540nm is directly proportional to glycerol concentration of the sample.



## ITEMS INCLUDED IN THE KIT

---

ITEM	DESCRIPTION	Cap Color	UNIT	QTY	STORAGE
Assay Plate, Plate A	96-well assay plate, blank	---	PLATE	2	-----
Dilution Buffer	12 ml	---	BOTTLE	1	4°C
Glycerol Reagent A	Reconstitute with 11.0 ml deionized water prior to use.	---	BOTTLE	1	4°C
Tray	For multi-channel pipetters, clear polyvinyl	CLEAR	EACH	2	-----
Glycerol standard	Glycerol @ 1mM [Dilute with 200 µl Dilution Buffer to make the 200 µM glycerol standard; see page 5 for recommended dilution scheme]	ORANGE	50 µl / VIAL	1	-20°C

### **Other equipment/reagents required but not provided with the kit:**


- Multi-channel Pipet , single channel pipet and pipet tips
- Plate reader with a filter of 540 nm
- Incubator at 37°C
- Large gauge needle
- Tubes for diluting standards

# ASSAY PROCEDURE

---

1. Prepare the glycerol standards as follows:

Briefly spin down the contents of the glycerol standard tube before reconstitution. Pipette 200  $\mu$ l of Dilution Buffer into the 1 mM glycerol standard tube provided and mix well by vortexing. This produces a diluted stock glycerol standard of 200  $\mu$ M. Pipette 125  $\mu$ l of dilution buffer into 6 tubes (not provided). Using the newly diluted stock glycerol solution, prepare a dilution series as depicted below. Mix each new dilution thoroughly before proceeding to the next. The 200  $\mu$ M stock dilution serves as the highest standard, and the dilution buffer serves as the zero standard.



2. Also at this time prepare the Glycerol Reagent A by adding 11.0 ml room temperature deionized water per bottle and gently invert. DO NOT VORTEX! Use a pipet to ensure that the powder is completely dissolved. Store at room temperature. If using a Reagent A solution previously prepared and stored at 2-8°C, also bring to room temperature. Make sure there is enough Reagent A from one solution to treat all the points in the assay. It may be necessary to combine solutions. Store in a light protected bottle. Reconstituted Glycerol Reagent A is stable for 60 days refrigerated (2-8°C); store any remaining solution refrigerated (2-8°C).
3. Add 20  $\mu$ l (or 10 - 25  $\mu$ l) of serum or plasma to a well of a blank plate for assessment of free glycerol. Add 30  $\mu$ l of dilution buffer to each well to total 50  $\mu$ l including serum or plasma sample. **THIS RESULTS IN A 2.5x DILUTION OF YOUR SAMPLE (20  $\mu$ l in 50  $\mu$ l).** Add 50  $\mu$ l of each standard to empty wells (use another plate, if necessary).
4. Add the reconstituted Glycerol Reagent A solution to one of the disposable trays provided in the kit. Add 50  $\mu$ l of Reagent A to each well of the glycerol plate. Gently, pipet up and down once to mix. Pop the bubbles using a large gauge needle or a clean pipet tip. The plate is then incubated at 25°C (room temperature) for 15 minutes.
5. The optical density of each well is then measured at 540 nm.

# GLYCEROL STANDARD CURVE

---


Generate standard curve: see example below

**DO NOT use this standard curve to generate your data. This is an example.]**

Subtract the OD value of the 0 $\mu$ M standard from all OD values including the standard curve.

Zero  
(blank) = .040

$\mu$ M Glycerol	OD	OD - blank
3.125	0.054	0.014
6.25	0.066	0.026
12.5	0.082	0.042
25	0.138	0.098
50	0.214	0.174
100	0.402	0.362
200	0.711	0.671



slope =	0.0034
intercept=	0.0075
$r^2$ =	0.9985

y = observed O.D. minus the blank

x = concentration of glycerol in  $\mu$ M

To calculate x for each y, (i.e. to change the observed O.D. into glycerol concentration) use the following equation:

y=(slope) times (x) plus intercept

$y=mx+b$  so  $x=(y-b)/m$

$x=(y - 0.0075)/0.003$  where 0.003= slope of the line and 0.0075= y intercept. Be careful to enter the proper sign for the y intercept value as it may be a negative number.

Any OD values greater than the highest standard (200  $\mu$ M), should be suspect. The compound should be re-assayed using a lower dose of the compound at treatment OR a dilute solution of the condition medium at the time of the assay.

**The  $R^2$  value should be equal or greater than 0.98 for the standard curve to be valid. Any  $R^2$  values below 0.98, must have the standard curve run again.**

Data are expressed as  $\mu$ M glycerol.

**REMEMBER TO ACCOUNT FOR THE DILUTION FACTOR IN STEP 3.**

APPENDIX A: PLATE LAYOUT

	A	B	C	D	E	F	G	H
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								

## APPENDIX B: PROCEDURE FLOWCHART

---

### Glycerol Detection

Add 20  $\mu$ l/well test sample and 30  $\mu$ l/well dilution buffer to one of the blank assay plates provided.  
Add 50  $\mu$ l/well diluted standard curve to empty wells.



Reconstitute Glycerol Reagent A.  
Add 50 $\mu$ l/well.



Incubate 15 minutes @ room temperature.



Measure the optical density of each well  
at 540 nm using a spectrophotometer  
plate reader.

**Reminder: Sample was diluted in Step 3**

