




# MOLECULAR ENDOCRINOLOGY

## MicroRNA 132 Regulates Nutritional Stress-Induced Chemokine Production through Repression of SirT1

Jay C. Strum, Jennifer H. Johnson, James Ward, Hongbo Xie, John Feild, Austin Hester, Alexander Alford and K. Michelle Waters

Mol. Endocrinol. 2009 23:1876-1884 originally published online Oct 9, 2009; , doi: 10.1210/me.2009-0117

To subscribe to *Molecular Endocrinology* or any of the other journals published by The Endocrine Society please go to: <http://mend.endojournals.org/subscriptions/>



## MicroRNA 132 Regulates Nutritional Stress-Induced Chemokine Production through Repression of SirT1

Jay C. Strum, Jennifer H. Johnson, James Ward, Hongbo Xie, John Feild, Austin Hester, Alexander Alford, and K. Michelle Waters

Departments of Metabolic Discovery Technology Group (J.C.S., J.H.J., H.X., A.H., A.A., K.M.W.), Computational Biology (J.W.), and Gene Cloning (J.F.), GlaxoSmithKline, Inc., Research Triangle Park, North Carolina 27709

Human adipose tissue secretes a number of proinflammatory mediators that may contribute to the pathophysiology of obesity-related disorders. Understanding the regulatory pathways that control their production is paramount to developing effective therapeutics to treat these diseases. Using primary human adipose-derived stem cells as a source of preadipocytes and *in vitro* differentiated adipocytes, we found IL-8 and monocyte chemoattractant protein-1 (MCP-1) are constitutively secreted by both cell types and induced in response to serum deprivation. MicroRNA profiling revealed the rapid induction of microRNA 132 (miR-132) in these cells when switched to serum-free medium. Furthermore, miR-132 overexpression was sufficient to induce nuclear factor- $\kappa$ B translocation, acetylation of p65, and production of IL-8 and MCP-1. Inhibitors of miR-132 decreased acetylated p65 and partially inhibited the production of IL-8 and MCP-1 induced by serum deprivation. MiR-132 was shown to inhibit silent information regulator 1 (SirT1) expression through a miR-132 binding site in the 3'-untranslated region of SirT1. Thus, in response to nutritional availability, induction of miR-132 decreases SirT1-mediated deacetylation of p65 leading to activation of nuclear factor- $\kappa$ B and transcription of IL-8 and MCP-1 in primary human preadipocytes and *in vitro* differentiated adipocytes. (*Molecular Endocrinology* 23: 1876–1884, 2009)

Historically, adipose tissue was thought to serve only as the storage site of excess energy; however, numerous studies have now identified a variety of proteins secreted from adipose, including both hormones and cytokines, which suggest a significant role for adipose in regulating whole-body energy metabolism and inflammation (1–9). These studies have led to the now generally accepted concept that adipose is a secretory organ with both endocrine and immunological activities (10). Obesity is a growing worldwide epidemic with type 2 diabetes and cardiovascular disease as common comorbidities. The development and progression of insulin resistance and atherosclerosis may result from a state of chronic inflammation within the adipose (11, 12). The pathogenesis of the inflammation is a matter of debate. Regardless of the cellular source of proinflammatory mediators within ad-

ipose, inhibition of their production may be therapeutically advantageous for the treatment of type 2 diabetes as well as cardiovascular disease. As a key transcription factor in inflammatory cascades in immune cells (13) and human adipocytes (14), nuclear factor- $\kappa$ B (NF $\kappa$ B) plays an important role in both innate and adaptive immunity. Chronic activation of NF $\kappa$ B *in vivo* has been shown to emulate the chronic inflammatory state observed with obesity and produces an insulin-resistant phenotype in IKK $\beta$  transgenic mice (15). However, the role of persistent activation of NF $\kappa$ B has not been directly assessed in human adipose. Furthermore, the mechanisms regulating the activation of NF $\kappa$ B and the production of proinflammatory mediators remain poorly characterized in human adipocytes.

MicroRNAs are small oligonucleotides that posttranscriptionally regulate gene expression by binding to the

ISSN Print 0888-8809 ISSN Online 1944-9917


Printed in U.S.A.

Copyright © 2009 by The Endocrine Society

doi: 10.1210/me.2009-0117 Received March 9, 2009. Accepted August 17, 2009.

First Published Online October 9, 2009

Abbreviations: CREB, cAMP response element binding protein; FBS, fetal bovine serum; IFN $\gamma$ , interferon- $\gamma$ ; MCP-1, monocyte chemoattractant protein-1; miR-132, microRNA 132; MSD, MesoScale Discovery; NF $\kappa$ B, nuclear factor- $\kappa$ B; PPAR $\gamma$ , peroxisome proliferator-activated receptor- $\gamma$ ; siRNA, small interfering RNA; SirT1, silent information regulator 1; 3'-UTR, 3'-untranslated region.



**FIG. 1.** Serum deprivation induces IL-8 and MCP-1. Human preadipocytes and *in vitro* differentiated adipocytes were incubated in serum-containing medium (■) or serum-free medium (◆) for the indicated times. After incubation, medium was collected, and the quantity of IL-8 and MCP-1 was determined using MSD assay kits. Proteins were quantitated by electrochemiluminescence detection and a MSD SECTOR Imager 6000. Data represent the mean of eight replicates ± SEM. These data are representative of three independent experiments.

3'-untranslated region (3'-UTR) of target gene sequences resulting in inhibition of transcription or translation (16, 17). Although over 700 microRNAs have been identified and various target predictions have been made, the function of these are just now being elucidated. Although the regulation of cell differentiation (18–20) by microRNAs and their dysregulation in various cancers (21, 22) is well documented, their role in signal transduction is less well understood.

To investigate the role of microRNAs in the production and secretion of proinflammatory mediators in human adipocytes, we profiled the expression of over 250 microRNAs in primary human preadipocytes and *in vitro* differentiated adipocytes under conditions of serum deprivation. We identified the induction of microRNA 132 (miR-132) and found it plays a key role in the response to serum deprivation in both preadipocytes and *in vitro* differentiated adipocytes. Furthermore, we demonstrate for the first time, the repression of silent information regulator 1 (SirT1) protein levels by miR-132 and propose a

pathway where miR-132 represses the SirT1-mediated deacetylation of p65 resulting in NFκB activation and IL-8 and monocyte chemoattractant protein-1 (MCP-1) production.


**Results**

**Serum deprivation induces IL-8 and MCP-1**

To assess the production of proinflammatory mediators by primary human preadipocytes and *in vitro* differentiated adipocytes, cells were cultured in the presence or absence of serum, and at various times, the medium was analyzed for the presence of interferon-γ (IFNγ), IL-1β, IL-10, IL-12p70, IL-13, IL-2, IL-4, IL-5, IL-8, TNFα, and MCP-1. Serum deprivation induced a significant increase in the chemokines IL-8 and MCP-1 in both cell types compared with cells grown in serum-containing medium (Fig. 1). No other cytokine tested was significantly produced from either cell type under these conditions (data not shown). These data show both cell types constitutively express IL-8 and MCP-1 when grown in serum-containing medium. However, the *in vitro* differentiated adipocytes secrete higher basal levels of IL-8 and MCP-1 and, in terms of chemokine production, are more sensitive to the stress of serum deprivation compared with the preadipocytes. These data suggest the regulation of chemokine production by nutrient availability.


**Serum deprivation induces miR-132**

To identify pathways involved in the induction of chemokines by serum deprivation, we investigated the effect of serum deprivation on the expression of microRNAs in human preadipocytes and *in vitro* differentiated human adipocytes. Initially, we measured the expression of 255 microRNAs using real-time PCR (data not shown). We identified an increase in miR-132 expression and, in subsequent studies, confirmed this observation in both cell types in response to serum deprivation (Fig. 2). Serum deprivation rapidly induced the biphasic expression of miR-132 with peaks at 1 and 24 h in both cell types.



**FIG. 2.** Serum deprivation induces miR-132 expression. *In vitro* differentiated human adipocytes and human preadipocytes were treated with serum-containing medium (■) or serum-free medium (◆) for the indicated times. Cells were lysed and total RNA extracted as described in *Materials and Methods*. Expression of miR-132 was determined using real-time PCR. Data represent the mean of four replicates and are representative of three independent experiments. Data are expressed as fold change relative to 0 h and were normalized to input of total RNA into the reverse transcription reaction.

The similarities between the miR-132 expression profile and chemokine production suggest a possible role for miR-132 in the regulation of chemokines in these cells.



**FIG. 3.** MiR-132 induces IL-8 and MCP-1. Human preadipocytes were plated ( $3 \times 10^4$  cells per well) and transfected with varying concentrations of pre-miR-132 or a scrambled oligonucleotide as a negative control (250 nM). After 48 h incubation, medium was collected and the quantity of IL-8 (white bars) and MCP-1 (black bars) was determined using MSD assay kits. Proteins were quantitated by electrochemiluminescence detection and a MSD SECTOR Imager 6000. Data are expressed as light units of cytokine and represent the mean of eight replicates  $\pm$  SEM. These data are representative of three independent experiments.

### miR-132 induces IL-8 and MCP-1

To gain insight into the possible relationship between induction of miR-132 and IL-8 and MCP-1 production, we explored the role of miR-132 by overexpressing miR-132 in preadipocytes and *in vitro* differentiated adipocytes. Transfection of pre-miR-132 in preadipocytes for 48 h induced a concentration-dependent increase of IL-8 and MCP-1 (Fig. 3). No other cytokine tested (IFN $\gamma$ , IL-1 $\beta$ , IL-10, IL-12p70, IL-13, IL-2, IL-4, IL-5, or TNF $\alpha$ ) was significantly produced from either cell type under these conditions (data not shown). The induction of both chemokines appeared biphasic with concentration of transfected pre-miR 132 with an initial plateau at 12.5–25 nM and subsequently at 200–250 nM. However, miR-132 was more efficacious at stimulating IL-8 production. We next investigated whether the ability of miR-132 to induce IL-8 and MCP-1 was unique to the pre-


dipocytes. Transfection of miR-132 for 48 h was similarly effective at inducing both chemokines in the *in vitro* differentiated adipocytes (data not shown). These data demonstrate that miR-132 overexpression is sufficient to induce IL-8 and MCP-1.

### Inhibitors of miR-132 partially inhibit induction of IL-8 and MCP-1 by serum deprivation


To determine whether blocking endogenous miR-132 affects the induction of IL-8 and MCP-1 by serum deprivation, human preadipocytes were transfected with inhibitors of miR-132. Inhibitors to miR-132 decreased the induction of IL-8 in a concentration-dependent manner (Fig. 4A). At 50 nM transfected inhibitor, there was an approximate decrease of 33% ( $P = 0.0002$ ). MiR-132 inhibitor was maximally effective at 250 nM (2-fold,  $P = 0.000002$ ) compared with serum withdrawal alone (Fig. 4A). Likewise, miR-132 inhibitors partially blocked the induction of MCP-1 by serum deprivation (Fig. 4B). The maximally effective concentration of miR-132 inhibitor was 100 nM with a decrease of 33% compared with cells grown in serum-free medium alone. This brought the level of MCP-1 down to about 2-fold greater than the constitutive production seen in cells grown in serum-containing medium. These changes were specific to the inhibitor used because the negative control had no significant effect on the production of IL-8 and MCP-1.

### Serum deprivation and miR-132 induces acetylation of p65 NF $\kappa$ B

NF $\kappa$ B plays a central role in the production of proinflammatory mediators, and IL-8 and MCP-1 are known transcriptional targets. Therefore, we initially assessed the role of NF $\kappa$ B by using 2-[(aminocarbonyl)amino]-5-(4-fluorophenyl)-3-thiophenecarboxamide (TPCA1), a potent and selective inhibitor of I $\kappa$ B kinase 2 (IKK-2). Pretreatment of either preadipocytes or *in vitro* differentiated adipocytes with 50 nM TPCA1 blocked the induction of IL-8 and MCP-1 by serum deprivation, suggesting NF $\kappa$ B activation is necessary for chemokine production (data not shown). We also measured translocation of the transcription factor from the cytosol to the nucleus in human preadipocytes. After 1 h incubation in serum-free medium, the translocation of NF $\kappa$ B from the cytosol into the nucleus was seen (data not shown). Likewise, after 24 h transfection, pre-miR-132 induced a similar response. To investigate the



**FIG. 4.** MiR-132 inhibitor partially blocks IL-8 and MCP-1 production. Human preadipocytes were transfected for 24 h in DMEM/F12 containing 10% FBS with varying concentrations of miR-132 inhibitor or a scrambled oligonucleotide negative control (500 nM). Cells were then incubated with fresh serum-containing medium or with serum-free medium for 24 h as indicated. After incubation, medium was collected, and the quantity of IL-8 (A) and MCP-1 (B) was determined using MSD assay kits. Proteins were quantitated by electrochemiluminescence detection and a MSD SECTOR Imager 6000. Data are expressed as light units of IL-8 and MCP-1 and represent the mean of eight replicates  $\pm$  SEM. These data are representative of three independent experiments.



**FIG. 5.** Serum deprivation and miR-132 overexpression induces acetylation of lysine 310 p65 and miR-132 inhibitors decrease acetylation. Human preadipocytes ( $1 \times 10^6$  cells) were plated in 25-cm<sup>2</sup> flasks. Cells were transfected for 24 h with a miR-132 inhibitor or scrambled oligonucleotide (A) or pre-miR-132 or a scrambled oligonucleotide control (B) before addition of serum-containing medium or serum-free medium for 1 h. After incubation, cells were lysed and protein was immunoprecipitated overnight at 4 C using anti-total p65. Proteins were resolved by SDS-PAGE and blotted with anti-acetyl lysine 310 antibody. Enhanced chemiluminescence detection was used. These data are representative of two independent experiments.

mechanism of NF $\kappa$ B activation, we measured the acetylation of lysine 310 of the p65 component of the complex. After 1 h in serum-free medium, an increase in endogenous acetylated lysine 310 p65 was seen in the preadipocytes (Fig. 5A) and *in vitro* differentiated adipocytes (data not shown) that was not seen in cells grown in serum-containing medium. Transfection of a miR-132 inhibitor decreased the acetylation of lysine 310 p65,


whereas a negative control oligonucleotide had no effect. Likewise, overexpression of miR-132 was sufficient to induce acetylation of lysine 310 p65 after 24 h transfection, whereas the transfection of a scrambled oligonucleotide control had no effect on acetylation of lysine 310 of p65 in cells grown in serum-free medium (Fig. 5B). Taken together, these data suggest the regulation of the acetylation/deacetylation pathway of NF $\kappa$ B by serum deprivation and miR-132.

### Serum deprivation and miR-132 regulate SirT1 expression

Because miR-132 appears to regulate the acetylation of p65, we searched for protein targets that are known to be involved in acetylation/deacetylation pathways. The nicotinamide adenine dinucleotide-dependent deacetylase SirT1 was found as a predicted target of miR-132 in MiRBase (23–25) (<http://microrna.sanger.ac.uk/>). Thus, we assessed the protein levels of SirT1 in human preadipocytes and *in vitro* differentiated adipocytes in response to serum deprivation and miR-132 overexpression. A significant decrease in SirT1 protein levels were observed at 1 and 3 h in cells incubated in serum-free medium compared with cells grown in serum-containing medium (Fig. 6A). Likewise, overexpression of miR-132 for 24 h decreased SirT1 protein levels relative to cells transfected with a scrambled oligonucleotide control (Fig. 6B) while transfection of antimirs to miR-132 blocked the decrease in SirT1 seen with serum deprivation (data not shown). The decrease in SirT1 protein with miR-132 overexpression suggests the direct translational repression of SirT1 by miR-132 and the possible regulation of SirT1-mediated deacetylation of p65 by miR-132.

### MiR-132 represses SirT1 expression through binding the 3'-UTR in the SirT1 gene

We first assessed whether serum deprivation or miR-132 overexpression affected SirT1 mRNA levels. Through gene expression analysis of SirT1 mRNA using real-time PCR, we found there was no change in SirT1 mRNA under any of these treatment conditions (data not shown). Thus, a luciferase reporter construct containing the 3'-UTR of the SirT1 gene and the predicted miR-132 binding site was prepared to assess whether the binding of miR-132 to the SirT1 3'-UTR mediates translational repression. After cotransfection of the reporter construct and pre-miR-132 or a scrambled oligonucleotide control in human preadipocytes, there was an approximately




**FIG. 6.** Serum deprivation and miR-132 overexpression decreases SirT1 protein levels. Human preadipocytes were plated in 12-well plates. Cells were incubated in serum-containing medium or serum-free medium for various times (A) or transfected with pre-miR-132 for 24 h (B). Cells were lysed, collected, and analyzed as described in *Materials and Methods*. The insets show the quantitation of SirT1 protein by densitometry. These data are representative of three independent experiments.

2-fold decrease in luciferase activity of the reporter construct in cells cotransfected with pre-miR-132 compared with the scrambled oligonucleotide control (Fig. 7). These data suggest miR-132 represses SirT1 translation through directly binding the response element in the SirT1 3'-UTR.

**Discussion**

Studies of microRNAs have yielded unique insights into the regulation of important biological processes.

Yamakuchi *et al.* (26) recently demonstrated the translational repression of SirT1 by miR-34a resulting in the activation of p53 and apoptosis, thus identifying miR-34a as a tumor suppressor. Our study is the first report to demonstrate the regulation of SirT1 by miR-132 and subsequent chemokine production through activation of NFκB. Previously, miR-132 was shown to be a transcriptional target of cAMP response element binding protein (CREB), which regulates neuronal morphogenesis through translational repression of p250GAP (27). Additionally, it's been



**FIG. 7.** miR-132 mediates translational repression of SirT1 through a binding site in SirT1 3'-UTR. Human preadipocytes ( $1 \times 10^5$  cells) were plated in 96-well plates and cotransfected with 200 ng pMIR-Report luciferase vector including the 3'-UTR of SirT1 containing the predicted miR-132 binding site and pre-miR-132 or scrambled oligonucleotide using an electroporator. Luciferase assays were performed using the dual-luciferase assay system 6–8 h after transfection. Data represent the mean of four biological replicates. The data are representative of two independent experiments.

shown to regulate the expression of methyl-CpG-binding protein (MeCP2) (28) and the angiotensin II type 1 receptor (AT1R) (29). MiR-132 was also reported to be induced by lipopolysaccharide in THP-1 cells (30), but a function was not ascribed. The biphasic expression pattern of miR-132 appears to parallel the biphasic activity profile of NF $\kappa$ B (31) and suggests miR-132 expression may be regulated by NF $\kappa$ B and represent a positive feedback loop. However, *in silico* promoter analysis failed to identify NF $\kappa$ B binding sites in the region of miR-132. As mentioned above, CREB has been shown to regulate miR-132 expression in neuronal cells (27). Additionally, CREB is regulated by high glucose in THP-1 cells where it has been shown to induce cyclooxygenase 2 (COX2), an enzyme that metabolizes arachidonic acid to prostaglandins (32). In our studies, CREB may regulate expression of miR-132 in response to the treatment conditions used where serum was absent and the concentration of glucose was relatively high (17.5 mM). Additional studies are needed to determine the exact mechanism of miR-132 transcription in human adipocytes.

SirT1 is highly regulated by nutrient availability. SirT1 protein levels *in vivo* are increased with starvation, fasting, and calorie restriction (33), whereas SirT1 protein decreases with age and senescence (33). Most interestingly, glucose appears to regulate SirT1 protein levels. Incubation of PC12 and HEK293 cells in the absence of both serum and glucose induces SirT1 protein expression through either an increase in transcription (34) or post-transcriptional regulation (35). In contrast, Nedachi *et al.* (36) showed low serum and high glucose represses SirT1 protein in C<sub>2</sub>C<sub>12</sub> cells. Likewise, Rodgers *et al.* (37) found high glucose mediated the posttranscriptional repression of SirT1 protein in mouse hepatocytes in conditions of low serum. Corroborating these findings, our study

shows SirT1 protein levels decrease when human preadipocytes and *in vitro* differentiated adipocytes are incubated in serum-free conditions and high glucose (17.5 mM) and extends these studies by proposing the posttranscriptional regulation of SirT1 by miR-132. Inhibition of SirT1 in response to high glucose in adipocytes should also result in activation of peroxisome proliferator-activated receptor- $\gamma$  (PPAR $\gamma$ ), which allows the cells to store glucose during times of excess energy (38). However, our data suggest nutrient availability may result in persistent activation of NF $\kappa$ B and chemokine production.

SirT1 has recently been implicated in the regulation of inflammation. Recently, Zhu *et al.* (39) showed resveratrol, a SirT1 activator, decreased TNF $\alpha$ -induced MCP-1 secretion in 3T3-L1 adipocytes. Pfluger *et al.* (40) showed overexpression of SirT1 in mice resulted in a lower recovery of IL-6 and TNF $\alpha$  in serum of transgenic mice fed a high-fat diet and an attenuated response to TNF $\alpha$ -induced NF $\kappa$ B activation in transgenic mouse embryonic fibroblasts. Thus, increased SirT1 activity appears to be antiinflammatory in mice. In contrast, when inhibited, SirT1 appears to be proinflammatory. Kwon *et al.* (41) showed HIV Tat binds SirT1 and inhibits the SirT1-mediated deacetylation of p65 NF $\kappa$ B resulting in hyperactivated NF $\kappa$ B, a greater transactivation efficiency and an exaggerated immune response. Likewise, Yang *et al.* (42) recently showed that cigarette smoke extract caused a dose- and time-dependent decrease in SirT1 activity and protein and a concomitant increase in NF $\kappa$ B-dependent proinflammatory mediator release in MonoMac 6 cells. More recently, studies using small interfering RNA (siRNA) to knock down SirT1 reported an increase in TNF $\alpha$ -induced MCP-1 and other proinflammatory genes in 3T3-L1 adipocytes (43). The authors went on to show that knockdown of SirT1 led to an increase in acetylation of lysine 310 of p65. Taken together with the findings reported here, these data demonstrate a decrease in SirT1 activity increases activation of NF $\kappa$ B and transcription of proinflammatory mediators.

Reversible acetylation of signaling proteins has emerged as a key posttranslational modification regulating activity of p53 and Forkhead proteins. Several reports have shown the reversible acetylation of NF $\kappa$ B (44–47). Furthermore, Yeung *et al.* (44) demonstrated the attenuation of NF $\kappa$ B activity by SirT1 and histone deacetylases leading to the enhancement of TNF $\alpha$ -induced apoptosis. Acetylation of lysine 310 of the p65 component of the NF $\kappa$ B complex has been shown to be necessary for the full activation of NF $\kappa$ B (45) and affect the duration and the selectivity of gene transcription (46). Ito (47) reported the acetylated lysine 310 p65 complex preferentially bound the NF $\kappa$ B binding site on the promoter region of IL-8 but not granulocyte-mac-

rophage colony-stimulating factor. We found miR-132 inhibitors blocked the accumulation of endogenous acetylated lysine 310 p65 and partially inhibited the induction of IL-8 and MCP-1 in serum-deprived cells. These data confirm the importance of acetylated lysine 310 p65 on the transcriptional activation of NF $\kappa$ B and extends it to conditions of nutrient availability in human preadipocytes and adipocytes. MiR-132 indirectly regulates acetylation of NF $\kappa$ B through repression of SirT1. SirT1 also binds to and regulates p300 (48), the acetyltransferase involved in NF $\kappa$ B activation. The decrease in SirT1 protein may have a possible dual effect on acetylation of NF $\kappa$ B through both inhibition of deacetylation and activation of acetylation.

Previously, IL-8 and MCP-1 were shown to be constitutively secreted by both human preadipocytes and isolated adipocytes implying a normal biological function for these chemokines (49–52). Excessive production through dysregulated NF $\kappa$ B is thought to play a pathogenic role in metabolic diseases. Indeed, circulating levels of IL-8 and MCP-1 are higher in obese than lean individuals and are associated with parameters of insulin resistance, atherosclerosis, and cardiovascular disease (53–59). Our data support the central role of NF $\kappa$ B in the inflammatory process in adipose and implies the miR-132 repression of SirT1-mediated deacetylation of p65 could be a mechanism of NF $\kappa$ B dysregulation contributing to the chronic inflammatory state in adipose tissue. Preventing the decrease in SirT1 protein should block the activation of NF $\kappa$ B and may be a therapeutic alternative to activating SirT1.

## Materials and Methods

### Cell culture and transfection

Primary human adipose-derived stem cells were collected from the abdominal sc adipose of a 37-yr-old female with a BMI of 23.3 (Zen-Bio Biological Systems and Services, Research Triangle Park, NC) and used as a source for preadipocytes and differentiated adipocytes. Cells (passages 4–8) were cultured in sc preadipocyte medium (Zen-Bio) until treatment. To differentiate into adipocytes, the sc preadipocytes were plated in preadipocyte medium and then after 24 h were differentiated in DMEM/F12 plus 10% fetal bovine serum (FBS) medium (Invitrogen, Carlsbad, CA) containing 200  $\mu$ M isobutylmethylxanthine, 20 nM dexamethasone, 20 nM GW347845X, a PPAR $\gamma$  agonist and 20 nM human insulin. After incubation for 7 d, the medium was changed to DMEM/F12 plus 10% FBS containing 20 nM insulin. All experiments on mature *in vitro* differentiated adipocytes were begun after d 14. Differentiation of cells was more than 90% as assessed by quantitating triglyceride accumulation using the Trinder assay (Sigma Chemical Co., St. Louis, MO) and real-time PCR analysis of fatty acid-binding protein 4 (FABP4 or ap2), PPAR $\gamma$ , and fatty acid synthase (FAS) mRNA, as markers of adipocyte differentiation. These results

showed less than 5% interwell variability among differentiated cultures of adipocytes. Serum deprivation was carried out by incubating cells in DMEM/F12 plus 1% BSA-fatty acid free. Experiments done to overexpress or inhibit miR-132 were done using Dharmafect 2 siRNA transfection reagent and miRIDIAN microRNA 132 mimic (pre-miR-132) UAACAGUCUACAGC-CAUGGUCG and mimic negative control UCACAACCUCUA-GAAAGAGUAGA (Dharmacon, LaFayette, CO) or miRCURY LNA knockdown probes CGACCATGGCTGTAGACTGTTA and scramble-miR control GTGTAACACGTCTATACGCCCA (Exiqon, Inc., Woburn, MA). The negative control oligonucleotides were tested at the highest concentration of the pre-miR-132 or inhibitor used in the experiment. Thus, we did not control for each concentration of microRNA used or transfect constant amounts of total microRNA. Uptake of microRNAs into cells was assessed using Cy-3-labeled pre-miR negative control and Cy-3-labeled anti-miR negative control (Ambion, Austin, TX) and Dharmafect 2 siRNA transfection reagent. Unless otherwise noted, all transfections were done in serum-containing medium.

### mRNA gene expression analysis

Total RNA was isolated from cells using Promega total RNA isolation kit (Promega, Madison, WI). RNA was reverse-transcribed using the ABI high-capacity cDNA archive kit (Applied Biosystems, Foster City, CA). Reactions containing 10 ng cDNA, SirT1 primers and probes (Assays on Demand; Applied Biosystems), and Universal PCR Master Mix (Applied Biosystems) were analyzed using a 7900HT sequence detector system (Applied Biosystems) as recommended by the vendor. The relative abundance of SirT1 mRNA was normalized to  $\beta$ -actin (Actb) mRNA.

### MicroRNA gene expression analysis

Total RNA was isolated by removing the treatment medium, washing the cells with PBS, and adding 1 ml TRIzol reagent per well. A phenol-chloroform extraction was performed, and 1.25 vol 100% ethanol was added to the aqueous phase. The microRNAs were isolated according to the manufacturer's protocol using a mirVana microRNA isolation kit from Ambion (Austin, TX). The microRNAs were then converted to cDNA, and gene expression analysis was done using the TaqMan microRNA assays (Applied Biosystems).

### NF $\kappa$ B translocation

Translocation of NF $\kappa$ B was assessed in human preadipocytes using the NF $\kappa$ B activation kit (ThermoScientific) and the Cello-mics ArrayScan II high-content imaging system. Briefly, human preadipocytes ( $3 \times 10^4$  cells per well) were plated into 96-well plates. After incubation overnight, cells were treated by incubating in serum-free medium or transfected with pre-miR-132 (100 nM) (Dharmacon) for 24 h. After treatment for various times, cells were fixed and stained with anti-NF $\kappa$ B antibody followed by an AlexaFluor 488-conjugated secondary antibody. The ratio of fluorescence in the cytosol and nucleus was calculated as a quantitative measure of NF $\kappa$ B translocation.

### Western blot analysis

Protein lysates were prepared using M-PER cell lysis buffer (Pierce, Rockford, IL) containing a protease cocktail inhibitor mix (Boehringer Mannheim, Mannheim, Germany). Protein lysates (10  $\mu$ g total protein) were resolved by electrophoresis on



10% or 4–12% SDS-polyacrylamide gels (Invitrogen, Carlsbad, CA) and transferred to nitrocellulose (Invitrogen). Proteins were detected by immunoblotting using anti-acetyl lysine 310 p65, anti-NF $\kappa$ B p65, anti-SirT1(C14H4), anti-GAPDH, antirabbit IgG horseradish peroxidase-linked antibody, anti-TFIIIB (Cell Signaling Technologies, Beverly, MA). All primary antibodies were used at 1:1000 dilution and incubated overnight at 4C. Antirabbit IgG HRP was used as the secondary antibody at 1:2000 dilution and detected using ECL reagents (Amersham, Piscataway, NJ).

### Plasmid construction

A 1.7-kb DNA fragment containing the SirT1 3'-UTR and the predicted miR-132 binding site was PCR amplified from human liver cDNA (QUICK-Clone cDNA; Clontech, Mountain View, CA). Primers were designed based on Genbank NM\_012238 and included a 5' *Sac*1 site and a 3' *Mlu*1 site to accommodate subcloning into the pMIR-REPORT plasmid (Ambion, Austin, TX).

### Luciferase assays

Primary cultures of human preadipocytes ( $1 \times 10^5$  cells per well) were plated in 96-well plates. After 24 h, 200 ng pMIR-Report luciferase vector including the 3'-UTR of SirT1, miR-132, or scrambled oligonucleotide were transfected using a Nucleofector (Amaxa, Walkersville, MD). Luciferase assays were performed using the dual-luciferase assay system (Promega) 6–8 h after transfection.

### Cytokine measurements

Media were collected at various times from preadipocytes and differentiated adipocytes to analyze secreted cytokines. MesoScale Discovery (MSD) assay kits were used according to the manufacturer's recommendations to quantitate IFN $\gamma$ , IL-1 $\beta$ , IL-10, IL-12p70, IL-13, IL-2, IL-4, IL-5, IL-8, TNF $\alpha$ , and MCP-1 (MesoScale Discovery, Gaithersburg, MD) protein levels using electrochemiluminescence detection and a MSD SECTOR Imager 6000.

### Acknowledgments

Address all correspondence and requests for reprints to: Jay C. Strum, 450 West Drive, Rm 21-232, Chapel Hill, North Carolina 27599-7295. E-mail: jaystrum@msn.com.

Disclosure Summary: All authors are employees of Glaxo-SmithKline (GSK). J.C.S., K.M.W., and J.F. own GSK stock options.

### References

- Schäffler A, Schölmerich J, Salzberger B 2007 Adipose tissue as an immunological organ: Toll-like receptors, C1q/TNFs and CTRPs. *Trends Immunol* 28:393–399
- Matarese G, La Cava A 2004 The intricate interface between immune system and metabolism. *Trends Immunol* 25:193–200
- Juge-Aubry CE, Henrichot E, Meier CA 2005 Adipose tissue: a regulator of inflammation. *Best Pract Res Clin Endocrinol Metab* 19:547–566
- Tilg H, Moschen AR 2006 Adipocytokines: mediators linking adipose tissue, inflammation and immunity. *Nat Rev Immunol* 6:772–783
- Schaffler A, Muller-Ladner U, Scholmerich J, Buchler C 2006 Role of adipose tissue as an inflammatory organ in human diseases. *Endocrinol Rev* 27:449–467
- Warne JP 2003 Tumor necrosis factor  $\alpha$ : a key regulator of adipose tissue mass. *J Endocrinol* 177:351–355
- Fain JN, Madan AK, Hiler ML, Cheema P, Bahouth SW 2004 Comparison of the release of adipokines by adipose tissue, adipose tissue matrix, and adipocytes from visceral and subcutaneous abdominal adipose tissues of obese humans. *Endocrinology* 145:2273–2282
- Hoch M, Eberle AN, Peterli R, Peters T, Sebock D, Keller U, Muller B, Linscheid P 2008 LPS induces interleukin-6 and interleukin-8 but not tumor necrosis factor- $\alpha$  in human adipocytes. *Cytokine* 41:29–37
- Chung S, Lapoint K, Martinez K, Kennedy A, Sandberg MB, McIntosh MK 2006 Preadipocytes mediate lipopolysaccharide-induced inflammation and insulin resistance in primary cultures of newly differentiated human adipocytes. *Endocrinology* 147:5340–5351
- Caspar-Bauguil S, Cousin B, Galinier A, Segafredo C, Nibbelink M, André M, Castella L, Pénicaud L 2005 Adipose tissues as an ancestral immune organ: site-specific change in obesity. *FEBS Lett* 579:3487–3492
- Wellen KE, Hotamisligil GS 2005 Inflammation, stress and diabetes. *J Clin Invest* 115:1111–1119
- Gerszten RE, Garcia-Zepeda EA, Lim YC, Yoshida M, Ding HA, Gimbrone Jr MA, Luster AD, Luscinskas FW, Rosenzweig A 1999 MCP-1 and IL-8 trigger firm adhesion of monocytes to vascular endothelium under flow conditions. *Nature* 398:718–723
- Caamaño J, Hunter CA 2002 NF- $\kappa$ B family of transcription factors: central regulators of innate and adaptive immune functions. *Clin Microbiol Rev* 15:414–429
- Creely SJ, McTernan PG, Kusminski CM, Fisher M, Da Silva NF, Khanolkar M, Evans M, Harte AL, Kumar S 2007 Lipopolysaccharide activates an innate immune system response in human adipose tissue in obesity and type 2 diabetes. *Am J Physiol Endocrinol Metab* 292:740–747
- Cai D, Yuan M, Frantz DF, Melendez PA, Hansen L, Lee J, Shoelson SE 2005 Local and systemic insulin resistance resulting from hepatic activation of IKK- $\beta$  and NF- $\kappa$ B. *Nat Med* 11:183–190
- Ambros V 2004 The function of animal microRNAs. *Nature* 431:350–355
- He L, Hannon GJ 2004 MicroRNAs: small RNAs with a big role in gene regulation. *Nat Rev Genet* 5:522–531
- Chen CZ, Li L, Lodish HF, Bartel DP 2004 MicroRNAs modulate hematopoietic lineage differentiation. *Science* 303:83–86
- Monticelli S, Ansel KM, Xiao C, Socci ND, Krichevsky AM, Thai TH, Rajewsky N, Marks DS, Sander C, Rajewsky K, Rao A, Kosik KS 2005 MicroRNA profiling of the murine hematopoietic system. *Genome Biol* 6:R71
- Esau C, Kang X, Peralta E, Hanson E, Marcusson EG, Ravichandran LV, Sun Y, Koo S, Perera RJ, Jain R, Dean NM, Freier SM, Bennett CF, Lollo B, Griffey R 2004 MicroRNA-143 regulates adipocyte differentiation. *J Biol Chem* 279:52361–52365
- Fabbri M, Garzon R, Andreeff M, Kantarjian HM, Garcia-Manero G, Calin GA 2008 MicroRNAs and noncoding RNAs in hematological malignancies: molecular, clinical and therapeutic implications. *Leukemia* 22:1095–1105
- Dalmay T 2008 MicroRNAs and cancer. *J Intern Med* 263:366–375
- Griffiths-Jones S, Saini HK, van Dongen S, Enright AJ 2008 miRBase: tools for microRNA genomics. *Nucleic Acid Res* 36(Database issue):D154–D158
- Griffiths-Jones S, Grocock RJ, van Dongen S, Bateman A, Enright AJ 2006 miRBase: microRNA sequences, targets and gene nomenclature. *Nucleic Acid Res* 34(Database issue):D140–D144
- Griffiths-Jones S 2004 The microRNA Registry. *Nucleic Acid Res* 32(Database issue):D109–D111
- Yamakuchi M, Ferlito M, Lowenstein CJ 2008 miR-34a repression of SIRT1 regulates apoptosis. *Proc Natl Acad Sci USA* 105:13421–13426

27. Vo N, Klein ME, Varlamova O, Keller DM, Yamamoto T, Goodman RH, Impey S 2005 A cAMP-response element binding protein-induced microRNA regulates neuronal morphogenesis. *Proc Natl Acad Sci USA* 102:16426–16431
28. Klein ME, Liyo DT, Ma L, Impey S, Mandel G, Goodman RH 2007 Homeostatic regulation of MeCP2 expression by a CREB-induced microRNA. *Nat Neurosci* 10:1513–1514
29. Elton TS 2008 MiR-132 regulates angiotensin II type 1 receptor expression through a protein coding region binding site. *Circulation* 118:S513
30. Taganov KD, Boldin MP, Chang KJ, Baltimore D 2006 NF- $\kappa$ B-dependent induction of microRNA miR-146, an inhibitor targeted to signaling proteins of innate immune responses. *Proc Natl Acad Sci USA* 103:12481–12486
31. Ladner KJ, Caligiuri MA, Guttridge DC 2003 Tumor necrosis factor-regulated biphasic activation of NF $\kappa$ B is required for cytokine-induced loss of skeletal muscle gene products. *J Biol Chem* 278:2294–2303
32. Shanmugam N, Gaw Gonzalo IT, Natarajan R 2004 Molecular mechanisms of high glucose-induced cyclooxygenase-2 expression in monocytes. *Diabetes* 53:795–802
33. Kwon HS, Ott M 2008 The ups and downs of SIRT1. *Trends Biol Sci* 609:1–9
34. Nemoto S, Fergusson MM, Finkel T 2004 Nutrient availability regulates SIRT1 through a forkhead-dependent pathway. *Science* 306:2105–2108
35. Kanfi Y, Peshti V, Gozlan YM, Rathaus M, Gil R, Cohen HY 2008 Regulation of SIRT1 protein levels by nutrient availability. *FEBS Lett* 582:2417–2423
36. Nedachi T, Kadotani A, Ariga M, Katagiri H, Kanzaki M 2008 Ambient glucose levels qualify the potency of insulin myogenic actions by regulating SIRT1 and FoxO3a in C<sub>2</sub>C<sub>12</sub> myocytes. *Am J Physiol Endocrinol Metab* 294:E668–E678
37. Rodgers JT, Lerin C, Haas W, Gygi SP, Spiegelman BM, Puigserver P 2005 Nutrient control of glucose homeostasis through a complex of PGC1- $\alpha$  and SIRT1. *Nature* 434:113–118
38. Picard F, Kurtev M, Chung N, Topark-Ngarm A, Senawong T, Machado De Oliveira R, Leid M, McBurney MW, Guarente L 2004 Sirt1 promotes fat mobilization in white adipocytes by repressing PPAR $\gamma$ . *Nature* 429:771–776
39. Zhu J, Yong W, Wu X, Yu Y, Lv J, Liu C, Mao X, Zhu Y, Xu K, Han X, Liu C 2008 Anti-inflammatory effect of resveratrol on TNF- $\alpha$ -induced MCP-1 expression in adipocytes. *Biochem Biophys Res Commun* 369:471–477
40. Pfluger PT, Herranz D, Velasco-Miguel S, Serrano M, Tschöp MH 2008 Sirt1 protects against high-fat diet-induced metabolic damage. *Proc Natl Acad Sci USA* 105:9793–9798
41. Kwon HS, Brent MM, Getachew R, Jayakumar P, Chen LF, Schnolzer M, McBurney MW, Marmorstein R, Greene WC, Ott M 2008 Human immunodeficiency virus type 1 Tat protein inhibits the SIRT1 deacetylase and induces T cell hyperactivation. *Cell Host Microbe* 3:158–167
42. Yang SR, Wright J, Bauter M, Seweryniak K, Kode A, Rahman I 2007 Sirtuin regulates cigarette smoke-induced proinflammatory mediator release via RelA/NF- $\kappa$ B in macrophages in vitro and in rat lungs in vivo: implications for chronic inflammation and aging. *Am J Physiol Lung Cell Mol Physiol* 292:L567–L576
43. Yoshizaki T, Milne JC, Imamura T, Schenk S, Sonoda N, Babendure JL, Lu JC, Smith JJ, Jirousek MR, Olefsky JM 2009 SIRT1 exerts anti-inflammatory effect and improves insulin sensitivity in adipocytes. *Mol Cell Biol* 29:1363–1374
44. Yeung F, Hoberg JE, Ramsey CS, Keller MD, Jones DR, Frye RA, Mayo MW 2004 Modulation of NF- $\kappa$ B-dependent transcription and cell survival by the SIRT1 deacetylase. *EMBO J* 23:2369–2380
45. Chen LF, Fischle W, Verdin E, Greene WC 2001 Duration of nuclear NF- $\kappa$ B action regulated by reversible acetylation. *Science* 293:1653–1657
46. Chen LF, Mu Y, Greene WC 2002 Acetylation of RelA at discrete sites regulates distinct nuclear functions of NF- $\kappa$ B. *EMBO J* 21:6539–6548
47. Ito K 2007 Impact of post-translational modifications of proteins on the inflammatory process. *Biochem Soc Trans* 35:281–283
48. Bouras T, Fu M, Sauve AA, Wang F, Quong AA, Perkins ND, Hay RT, Gu W, Pestell RG 2005 SIRT1 deacetylation and repression of p300 involves lysine residues 1020/1024 within the cell cycle regulatory domain 1. *J Biol Chem* 280:10264–10276
49. Bruun JM, Pedersen SB, Richelsen B 2001 Regulation of interleukin 8 production and gene expression in human adipose tissue *in vitro*. *J Clin Endocrinol Metab* 86:1267–1273
50. Gerhardt CC, Romero IA, Canello R, Camoin L, Strosberg AD 2001 Chemokines control fat accumulation and leptin secretion by cultured human adipocytes. *Mol Cell Endocrinol* 175:81–92
51. Bruun JM, Lihn AS, Madan AK, Pedersen SB, Schiøtt KM, Fain JN, Richelsen B 2004 Higher production of IL-8 in visceral vs. subcutaneous adipose tissue. Implication of nonadipose cells in adipose tissue. *Am J Physiol Endocrinol Metab* 286:E8–E13
52. Christiansen T, Richelsen B, Bruun JM 2005 Monocyte chemoattractant protein 1 is produced in isolated adipocytes, associated with adiposity and reduced after weight loss in morbid obese subjects. *Int J Obes* 29:146–150
53. Wang B, Jenkins JR, Trayhurn P 2005 Expression and secretion of inflammation-related adipokines by human adipocytes differentiated in culture: integrated response to TNF- $\alpha$ . *Am J Physiol Endocrinol Metab* 288:E731–E740
54. Charo IF, Taubman MB 2004 Chemokines in the pathogenesis of vascular disease. *Circ Res* 95:858–866
55. Sartipy P, Loskutoff DJ 2003 Monocyte chemoattractant protein 1 in obesity and insulin resistance. *Proc Natl Acad Sci USA* 100:7265–7270
56. Bruun JM, Lihn AS, Pedersen SB, Richelsen B 2005 Monocyte chemoattractant protein-1 release is higher in visceral than subcutaneous human adipose tissue (AT): implication of macrophages resident in the AT. *J Clin Endocrinol Metab* 90:2282–2289
57. Moreau M, Brocheriou I, Petit L, Ninio E, Chapman MJ, Rouis M 1999 Interleukin-8 mediates downregulation of tissue inhibitor of metalloproteinase-1 expression in cholesterol-loaded human macrophages: relevance to stability of atherosclerotic plaque. *Circulation* 99:420–426
58. Straczkowski M, Dzienis-Straczkowska S, Stępień A, Kowalska I, Szelachowska M, Kinalska I 2002 Plasma interleukin-8 concentrations are increased in obese subjects and related to fat mass and tumor necrosis factor- $\alpha$  system. *J Clin Endocrinol Metab* 87:4602–4606
59. Zozulińska D, Majchrzak A, Sobieska M, Wiktorowicz K, Wierusz- Wysocka B 1999 Serum interleukin-8 level is increased in diabetic patients. *Diabetologia* 42:117–118

